



# Transport Overlay



Transport overlays are available in many different versions and customised to the respective transport stretcher. Depending on the model, there are versions equipped with carrying straps, which serve to transfer the patient including overlay to another base frame, or side bolsters for improved positioning of the patient. All versions have one thing in common: The PUNL® cover and the lining with viscoelastic foam (see detailed description below).

## MOLTOSLOW



The viscoelastic thermal foam reacts to pressure and warmth and optimises the adaption to the body, thus creating a full-surface body support automatically adjusting to the body contours in order to ensure optimum pressure distribution. The shear stress is significantly reduced, whilst the sensory pressure relief additionally ensures some ‚fixation‘ of the patient. The risk of suffering pressure sores is reduced to a minimum. Example of a transport overlay

Example

