

## Product- & Pricelist MOMOSAN for Veterinary Wound Care

<b>Item no. 61123</b>		
10 pc.	MOMOSAN® white wound dressing, large sheets, low-germ,	59x49x2cm
<b>Item no. 61121</b>		
20 pc.	MOMOSAN® white wound dressing, large sheets, low-germ,	59x49x1cm
<b>Item no. 61601</b>		
143 pc.	MOMOSAN® white wound dressing, small sheets, low-germ,	15x10x2cm
<b>Item no. 61602</b>		
286 pc.	MOMOSAN® white wound dressing, small sheets, low-germ,	15x10x1cm
<b>Item no. 61605</b>		
30 pc.	MOMOSAN® white wound dressing, small sheets, <b>sterile</b> ,	15x10x0,5cm
<b>Item no. 61901</b>		
20 pc.	MOMOSAN® white wound dressing, small sheets, <b>sterile</b> ,	15x10x1cm
<b>Item no. 61903</b>		
10 pc.	MOMOSAN® white wound dressing, small sheets, <b>sterile</b> ,	24x16x1cm
<b>Item no. 61905</b>		
14 pc.	MOMOSAN® white tamponade, <b>sterile</b> ,	200x2,5x0,4cm
<b>Item no. 61132</b>		
12 pc.	MOMOSAN® white bandages, low-germ,	300x5x0,3cm
<b>Item no. 61128</b>		
6 pc.	MOMOSAN® white bandages, low-germ,	300x10x0,3cm
<b>Item no. 61666</b>		
10 pc.	MOMOSAN® black wound pad, climatize, low-germ,	55x45x2cm
<b>Item no. 61667</b>		
2 pc.	MOMOSAN® black wound pad, climatize, low-germ,	24x16x5cm
<b>Item no. 61920</b>		
30 pc.	MOMOSAN® white <b>Woundcleaner</b> , <b>sterile</b> ,	5x5x2cm
<b>Item no. 61921</b>		
30 pc.	MOMOSAN® black <b>Woundcleaner</b> , <b>sterile</b> ,	5x5x2cm
<b>Item no. 61922</b>		
15 pc.	MOMOSAN® white <b>Woundcleaner</b> , <b>sterile</b> ,	5x5x2cm
15 pc.	MOMOSAN® black <b>Woundcleaner</b> , <b>sterile</b> ,	5x5x2cm
<b>Item no. 61520 Set VET Large Animals</b>		
10 pc.	MOMOSAN® white wound dressing, low-germ,	59x49x2cm
10 pc.	MOMOSAN® black wound pad, climatize,	55x45x2cm
14 pc.	MOMOSAN® white tamponade, <b>sterile</b> ,	200x2,5x0,4cm
12 pc.	MOMOSAN® white tamponade, low-germ,	300x10x0,3cm
20 pc.	MOMOSAN® white wound dressing, <b>sterile</b> ,	24x16x1cm
30 pc.	MOMOSAN® black Woundcleaner, <b>sterile</b> ,	5x5x2cm
26 pc.	MOMOSAN® white wound dressing, low-germ	15x10x2cm
<b>Item no. 61521 Set VET Small Animals</b>		
20 pc.	MOMOSAN® white wound dressing, low-germ,	59x49x1cm
10 pc.	MOMOSAN® black wound pad, climatize,	55x45x2cm
14 pc.	MOMOSAN® white tamponade, <b>sterile</b> ,	200x2,5x0,4cm
24 pc.	MOMOSAN® white tamponade, low-germ,	300x5x0,3cm
20 pc.	MOMOSAN® white wound dressing, <b>sterile</b> ,	15x10x1cm
30 pc.	MOMOSAN® white wound dressing, <b>sterile</b> ,	15x10x0,5cm
30 pc.	MOMOSAN® white Woundcleaner, <b>sterile</b> ,	5x5x2cm
52 pc.	MOMOSAN® white wound dressing, low-germ	15x10x1cm
<b>Item no. 61500 Set VET Standard</b>		
5 pc.	MOMOSAN® white wound dressing, large sheets, low-germ,	59x49x2cm
5 pc.	MOMOSAN® white wound dressing, large sheets, low-germ,	59x49x1cm
3 pc.	MOMOSAN® white tamponade, low-germ,	300x5x0,3cm
<b>Item no. 61580 Set VET Sheets „XXL“</b>		
3 pc.	MOMOSAN® white wound dressing, large sheets, low-germ,	200x50x2cm
<b>Item no. 61570 Set VET climatized sheets</b>		
6 pc.	MOMOSAN® black	50x45x2cm
6 pc.	MOMOSAN® black	50x45x1cm

More different sizes and contents of package can be found in the MOMOSAN® productlist for human medicine: [www.momosan.com/e\\_productlist.pdf](http://www.momosan.com/e_productlist.pdf)  
Application referred to instruction for use.  
The prices indicated are exclusive VAT.



Simple and effective  
Wound Care for

Large animals  
Small animals  
Wild animals

Producer:



**MOLTOPLAST®** healthcare

A-6020 Innsbruck · Egger-Lienz-Straße 130  
phone. +43 (0)512 / 58 54 07 · Fax: +43 (0)512 / 58 54 07-18  
med@moltoplast.com · www.momosan.com



## Product Information

### MOMOSAN® – veterinary wound care

...is a polyurethane wound dressing in honeycomb structure with controlled absorbing and pressure-relieving effect

...effective without using medicine and active substances

...causes controlled exudate absorption whilst maintaining a moist wound bed

...removes and absorbs dead tissue

...prevents new contamination

...promotes the formation of granulation tissue and epithelialisation upward and inside


...does not dissolve inside the wound

...no creams, lotions or ointments required

...can be cut to size and formed as desired

...inexpensive and easy to handle


...both sterile and non-sterile versions available\*



MOMOSAN white



MOMOSAN black



\* The medical practitioner has to decide about the respective use of sterile or non-sterile products.



## Application Areas


### MOMOSAN® white – veterinary wound care

- All kinds of wounds in horses, cattle, sheep, pigs, dogs, cats and wildlife
- Superficial wounds, deep wounds, pocket wounds
- Wounds with small openings/fistulas
- Traumatic wounds
- Chronic wounds
- Infected wounds
- Décollement (detachment of the skin)
- Abscesses
- Burns
- Stabbing and cutting injuries
- Wound dehiscence
- Covering of fresh stitches
- Talus, ankle joint and ankle injuries
- Prevention of uncontrolled tissue growth in case of injuries

### MOMOSAN® black – veterinary wound care

- Treatment of open, septic/deeply suppurating wounds with high exudate management
- Débridement
- Dressing for vacuum therapy
- climatized base upholstery
- cage padding

### Wound care - example:



Information of the correct application can be found in the corresponding instruction for use of MOMOSAN® products.

### Suture insufficiency after pan-arthrodesis - dog



### Injury - horse



### Surgical suture - dog



### Injury - horse



Result after 6 days MOMOSAN white

# Success in animal healthcare